

Clamp Electrodes

- Universal connection is suitable for use with 3.2 mm straight pin, 4.0 mm banana or snap lead termination leadwires.
- All of the metallic surfaces are Silver–Silver Chloride (Ag/AgCl) plated to provide maximum signal quality.
- Color coded IEC or AHA protocols.
- The lengths of clamp electrodes are 140 mm for adult size, 89 mm for pediatric size.
- Supplied in package of 4 pieces set.

Adult	IEC	CEAI1	Package 4 pcs./set, Red, Green, Yellow, Black
Adult	AHA	CEAA1	Package 4 pcs./set, Red, Green, White, Black
Pediatric	IEC	CEAI2	Package 4 pcs./set, Red, Green, Yellow, Black
Pediatric	AHA	CEAA2	Package 4 pcs./set, Red, Green, White, Black



Suction Chest Electrode

- Suction chest electrodes with universal connection, are suitable for use with 3.2 mm straight pin, 4.0 mm banana or snap lead termination leadwires.
- All of the metallic surfaces are Silver–Silver Chloride (Ag/AgCl) plated to provide maximum signal quality.
- The outer diameter of adult suction bowl is 24.0 mm. The diameter of the silicone bulb is 30.0 mm. The diameter of pediatric suction bowl is 17.0 mm. The diameter of the silicone bulb is 21.4 mm.
- Supplied in package of 6 pieces set.

Adult	SE1	Package 6 pcs./set
Pediatric	SE2	Package 6 pcs./set

Banana to Snap Adaptor

- Banana to Snap Adaptors allows you to easily connect 3.2 mm straight pin and 4.0 mm banana plug leadwires to disposable ECG electrodes with 4.0 mm snap studs.
- Banana to Snap Adaptors supplied in package of 10 pieces set.



- BS3 3.2 mm straight pin to snap adaptor, Grey, Package 10 pcs./set
- BS4 4.0 mm banana to snap adaptor, White, Package 10 pcs./set

Banana to Alligator Clip Adaptor

- Allows you to easily connect 3.0 mm or 4.0 mm banana plug leadwires to disposable ECG electrodes with tab connectors.
- Supplied in package of 10 pieces set.

BA3	3.0 mm banana to alligator adaptor, 10 pcs./set
BA4	4.0 mm banana to alligator adaptor, 10 pcs./set



Banana to Grabber Adaptor

- Converts 3.2 mm straight pin and 4.0 mm banana plug leadwires to disposable ECG electrodes with 4.0 mm snap studs.
- Supplied in package of 10 pieces set.

BG3I/10	3.2 mm straight pin to grabber adaptor, Grey, IEC color code, Package	10 pcs./set
BG3A/10	3.2 mm straight pin to grabber adaptor, Grey, AHA color code, Package	10 pcs./set
BG4/10	4.0 mm banana to grabber adaptor, White, Package	10 pcs./set

