


Diagnostic EKG Cables & Accessories

Creating an Order Code



- 00 FMT E100 Monolithic (non detachable leads)
- 01 FMT E101 All detachable leads
- 03 FMT E103 5+5 detachable leads (GE Multi Link)
- 04 FMT E104 5+5 detachable leads (GE Multi Link 12SL)
- 05 FMT E105 5+5 detachable leads (Philips)
- 06 FMT E106 4+6 detachable leads (Medtronic Physio Control)
- 07 FMT E107 4 Fixed 6 Replaceable Leads (Schiller)

- B 4.0 mm isolated banana
- P 3.2 mm straight pin
- S Snap
- G Grabber
- 0 Trunk cable without leadwire

To determine the part code enter specified parameter values instead of X.

- 1 FMT Diagnostic EKG cables are monolithic or with replaceable leads compatible with various trunk cable types for different EKG devices. Available options may be determined from "Yoke Types" section.

Full part number is valid for trunk cable and grouped leads. If you would like to buy the trunk cable without leads, just write "0" value for patient end -field 4-, after filling other sections.

- 2 You can find cables for various brands and models, resistor values in "Compatible Models" section.

- 3 4 Patient end connections can be selected from banana, straight pin, grabber or snap options. You can choose Patient End code and IEC or AHA color code from "Leadwires Patient Ends" section.

EKG Cables have 1.80 m trunk cable length (Ø6.0 mm), 80 cm chest and 100 cm extremity lead lengths (Ø3.0 mm). Other lengths are available upon request.

Yoke Types

E100 type yoke



One piece 10 lead (non detachable leads).

E101 type yoke



10 lead (all detachable leads)

E103 type yoke



10 lead (5+5 detachable leads),
GE MULTI-LINK type yoke.

E104 type yoke



10 lead (5+5 detachable leads),
GE MULTI-LINK 12SL type yoke.

E105 type yoke



10 lead (5+5 detachable leads),
PHILIPS type yoke.

E106 type yoke



10 lead (4+6 detachable leads),
Medtronic/Physio Control type yoke.


E107 type yoke



10 lead (4+6 detachable leads),
Schiller type yoke.

Diagnostic EKG Cables & Accessories

Compatible Brands & Models



15 pin D-sub connector with screws




E100/DNY/Z
E101/DNY/Z

E100/DRY/Z
E101/DRY/Z

10 kΩ

Philips/HP M1770A(Pagewriter 200i-300i), M1771A (Pagewriter 200), M1772A (Pagewriter 100), Cable OEM M3703C. **Bionet** All Models. **BTL** All Models. **Aspel** All Models. **Innomed** All Models. **Mortara** ELI150, ELI210, ELI250, ELI350.




15 pin D-sub connector with screws



E100/SRY/Z
E101/SRY/Z

10 kΩ

Schiller AT, AT1, AT2, AT3, AT5, AT6, AT10, AT60, A104 PC, CS100, CS6, CS6/12, AT101, AT102, AT104. **Esaote** P80, P120. **Kontron** P80.



15 pin D-sub connector with screws



E100/NRY/Z
E101/NRY/Z

E100/NRY/Z
E101/NRY/Z

10 kΩ

Nihon Kohden Cardiofax C6151, 6353, 8270, 8350, 8370, 8420, 8423, 8450, 8453, 8820, 8823, 8830, 9620, 9010, 9020, 9022, 9110, 9130, 9132. **Dong Jiang** All Models. **Biocare** All Models. **Cardioline** AR600, AR1200, ELAN, ELAN1100, Cardiette Excel 103, 106, 112, 212, Deadalus Models, Microruler & Autoruler Models.




15 pin D-sub connector



E100/GRY/Z
E101/GRY/Z
E103/GRY/Z

10 kΩ

GE Microsmart, MAC500 MAC1000, MAC1100, MAC1200ST, MAC400, MAC600, MAC800, MAC1600.



15 pin D-sub connector with screws



E100/ERY/Z
E101/ERY/Z

10 kΩ

Edan SE-1200 Express, SE-12, SE-12 Express, SE-601, SE-3, SE-1. **Mindray** DECG-03A.



16 pin circular connector



E100/CRY/Z
E101/CRY/Z

4.7 kΩ

Cardioline Delta 1, 3, 30, 60, Delta 1 plus, Delta 3 plus, Delta 60 plus. **Suzuken Kenz** 103, 106.



15 pin D-sub connector with latch on the top



E100/KRY/Z
E101/KRY/Z

10 kΩ

Fukuda ME KP-5000.



15 pin D-sub connector with 45 degree angle and locking plastic latches



E100/HRY/Z
E101/HRY/Z

10 kΩ

Bosch 603D, 603ER, 603S, 606. **Customed**, **Hellige** EK36, EK41, EK43, EK53, EK56, EK403, 413, 512, E512P, Cardiosmart, Cardiosys. **Siemens** Sicard 440/440S, 460/460S, 740, Sicard 803, Sirecust 701, Cathcor, Cardiostat 1,3T, 703, 710, Cardinex 3T, 4, 6T, 703, 710, Megacard Elema A13, R, R/E, Miracard1. **Hörmann** Bioset9000.



15 pin D-sub connector with locking metal latches



E100/FNY/Z
E101/FNY/Z

-

E100/FRY/Z
E101/FRY/Z

10 kΩ

Fukuda Denshi Autocardiner, Cardiomax, FX2111, FX7411, FX7402, FX7202, FX2155, FX3010, FX7102.