


Diagnostic EKG Cables & Accessories

Compatible Brands & Models



15 pin D-sub connector with screws




E100/DNY/Z
E101/DNY/Z

E100/DRY/Z
E101/DRY/Z

10 kΩ

Philips/HP M1770A(Pagewriter 200i-300i), M1771A (Pagewriter 200), M1772A (Pagewriter 100), Cable OEM M3703C. **Bionet** All Models. **BTL** All Models. **Aspel** All Models. **Innomed** All Models. **Mortara** ELI150, ELI210, ELI250, ELI350.




15 pin D-sub connector with screws



E100/SRY/Z
E101/SRY/Z

10 kΩ

Schiller AT, AT1, AT2, AT3, AT5, AT6, AT10, AT60, A104 PC, CS100, CS6, CS6/12, AT101, AT102, AT104. **Esaote** P80, P120. **Kontron** P80.



15 pin D-sub connector with screws



E100/NNY/Z
E101/NNY/Z

E100/NRY/Z
E101/NRY/Z

10 kΩ

Nihon Kohden Cardiofax C6151, 6353, 8270, 8350, 8370, 8420, 8423, 8450, 8453, 8820, 8823, 8830, 9620, 9010, 9020, 9022, 9110, 9130, 9132. **Dong Jiang** All Models. **Biocare** All Models. **Cardioline** AR600, AR1200, ELAN, ELAN1100, Cardiette Excel 103, 106, 112, 212, Deadalus Models, Microruler & Autoruler Models.



15 pin D-sub connector



E100/GRY/Z
E101/GRY/Z
E103/GRY/Z

10 kΩ

GE Microsmart, MAC500 MAC1000, MAC1100, MAC1200ST, MAC400, MAC600, MAC800, MAC1600.



15 pin D-sub connector with screws



E100/ERY/Z
E101/ERY/Z

10 kΩ

Edan SE-1200 Express, SE-12, SE-12 Express, SE-601, SE-3, SE-1. **Mindray** DECG-03A.



16 pin circular connector



E100/CRY/Z
E101/CRY/Z

4.7 kΩ

Cardioline Delta 1, 3, 30, 60, Delta 1 plus, Delta 3 plus, Delta 60 plus. **Suzuken Kenz** 103, 106.



15 pin D-sub connector with latch on the top



E100/KRY/Z
E101/KRY/Z

10 kΩ

Fukuda ME KP-5000.



15 pin D-sub connector with 45 degree angle and locking plastic latches



E100/HRY/Z
E101/HRY/Z

10 kΩ

Bosch 603D, 603ER, 603S, 606. **Customed**, **Hellige** EK36, EK41, EK43, EK53, EK56, EK403, 413, 512, E512P, Cardiosmart, Cardiosys. **Siemens** Sicard 440/440S, 460/460S, 740, Sicard 803, Sirecust 701, Cathcor, Cardiostat 1,3T, 703, 710, Cardinex 3T, 4, 6T, 703, 710, Megacard Elema A13, R, R/E, Miracard1. **Hörmann** Bioset9000.



15 pin D-sub connector with locking metal latches



E100/FNY/Z
E101/FNY/Z

-

E100/FRY/Z
E101/FRY/Z

10 kΩ

Fukuda Denshi Autocardiner, Cardiomax, FX2111, FX7411, FX7402, FX7202, FX2155, FX3010, FX7102.