

Probes with hose connection



Probe type **AFNOR NF-S 90-116 EASYFIX®** for compressed medical gases with hose connector Ø 6 mm and for Vacuum with hose connector Ø 7.5 mm.




Probe type **UNI 9507** for compressed medical gases with hose connector Ø 6 mm and for Vacuum with hose connector Ø 7.5 mm.



Probe type **NIST EN ISO 5359** for compressed medical gases with hose connector Ø 6 mm and for Vacuum with hose connector Ø 7.5 mm.



Probe type **BS 5682** for compressed medical gases with hose connector Ø 6 mm and for Vacuum with hose connector Ø 7.5 mm.



Probe type **DIN 13260-2** for compressed medical gases with hose connector Ø 6 mm and for Vacuum with hose connector Ø 7.5 mm.



Probe type **SS 875 24 30** for compressed medical gases with hose connector Ø 6 mm and for Vacuum with hose connector Ø 7.5 mm.



Probe type **JIS T 7101** for compressed medical gases with hose connector Ø 6 mm and for Vacuum with hose connector Ø 7.5 mm.



Probe type **SANS 1409** for compressed medical gases with hose connector Ø 6 mm and for Vacuum with hose connector Ø 7.5 mm.



Probe type **OHMEDA** for compressed medical gases with hose connector Ø 6 mm and for Vacuum with hose connector Ø 7.5 mm.